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L14

## 2.1.14 Anion Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Nucleosil SB (strong anion exchanger)

The anion exchanger Nucleosil SB is used for ion exchange chromatography. It is a strongly basic anion exchanger with a quaternary ammonium modification and a capacity of about 1 meq/g.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... <b>E150NS050</b>	Nucleosil 100-5-SB	SAX	5 µm	spherical	100 Å	350 m <sup>2</sup> /g	10	no	MACHEREY-NAGEL
.... <b>E150NS100</b>	Nucleosil 100-10-SB	SAX	10 µm	spherical	100 Å	350 m <sup>2</sup> /g	10	no	MACHEREY-NAGEL

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## 2.1.14 Anion Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Separation of RNA derivatives by anion exchange

#### Nucleosil 100-10 SB

Dimension: 300 x 4 mm  
 Eluent: 0.2M ammonium acetate pH 7.5  
 Flow rate: 1 ml/min  
 Pressure: 12400 - 13800 kPA  
 Detection: UV, 254 nm  
 Temperature: ambient  
 Injection: 3 µl


Author(s): Alonso J.  
 Source: J. Liquid Chromatogr. 8 (1985) 299 - 313

Substances: cytidine, RNA  
 Fig.: Elution profile of a mixture of cytidine (Cyd), cytidine 2'- and 3'-phosphate (2'-CMP, 3'-CMP) and cytidine 2':3'-phosphate (2':3'-CMP).

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*The chemicals mentioned in the application notes are neither manufactured nor sold by AfC or MN. The applications described in literature or in our data bases are just for your information. Therefore we can not guarantee for any test or correctness of the applications.*

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L14

## 2.1.14 Anion Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Partisil SAX

Partisil was one of the first commercial available type A silica packings on the HPLC market. The benefit of these packings compared to others was the good batch to batch reproducibility of those stationary phases. It is an irregular shaped support. Several different bonding types are existing. There are three different C18 packings available. Whereas Partisil ODS is a low density monomeric C18 phase with no endcapping, the Partisil ODS(2) is polymeric bonded and Partisil ODS(3) has an additional endcapping combined with a polymeric bonding. The Partisil PAC phase is an alkyl bonded phase with a mixed Amino/Cyano endcapping.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... <b>D650PA100</b>	Partisil SAX	SAX	10 µm	irregular	85Å	350 m <sup>2</sup> /g	10	no	Whatman

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L14

## 2.1.14 Anion Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Zorbax SAX

The classical Zorbax product line is typical for family of spherical first generation (type A) packings. The stationary phases show a high hydrophobicity and the support has an excellent mechanical strength. AfC offers the complete product line of the classical Zorbax supports in all column dimensions.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... <b>D150ZX050</b>	Zorbax SAX	SAX	5 µm	spherical	70 Å	300 m <sup>2</sup> /g	-	no	Agilent Technologies

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L14

## 2.1.14 Anion Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Sphere-Image 80 SAX

The technical specifications of the Sphere-Image 80 series of packing materials are nearly equivalent to those of Spherisorb 80, a high-performance packing material which has existed for decades and is very well-known. The Sphere-Image 80 packing material is also available with all of the typical modifications. The cation and anion exchangers Sphere-Image SCX and SAX are also available. The most popular materials in the Sphere-Image series are ODS 1 (not endcapped) and end-capped ODS 2. All of these packing materials exhibit excellent batch-to-batch reproducibility.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... <b>D150SI050</b>	Sphere-Image SAX	SAX	5 µm	spherical	80 Å	200 m <sup>2</sup> /g	-	no	Phenomenex
.... <b>D150SI100</b>	Sphere-Image SAX	SAX	10 µm	spherical	80 Å	200 m <sup>2</sup> /g	-	no	Phenomenex

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• Bulk material available in 10 g and 100 g quantities.

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## 2.1.14 Cation Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Nucleosil SA (strong cation exchanger)

The cation exchanger Nucleosil SA is used for ion exchange chromatography. It is a strongly acidic cation exchanger with sulphonic acid modification and a capacity of about 1 meq/g.

- Bulk material available in 10 g and 100 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.

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## 2.1.14 Cation Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Separation of nucleosides by HPLC


#### Nucleosil 100-10 SA

Dimension: 300 x 4 mm  
Eluent: 0.2M (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> pH 4.0, adjusted with HCl  
Flow rate: 0.25 cm/sec  
Pressure: 130 bar  
Detection: UV, 254 nm  
Temperature: ambient  
Injection: 3 µl

Author(s): Macherey-Nagel  
Source: LC department

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L9

## 2.1.14 Cation Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Partisil SCX

Partisil was one of the first commercial available type A silica packings on the HPLC market. The benefit of these packings compared to others was the good batch to batch reproducibility of those stationary phases. It is an irregular shaped support. Several different bonding types are existing. There are three different C18 packings available. Whereas Partisil ODS is a low density monomeric C18 phase with no endcapping, the Partisil ODS(2) is polymeric bonded and Partisil ODS(3) has an additional endcapping combined with a polymeric bonding. The Partisil PAC phase is an alkyl bonded phase with a mixed Amino/Cyano endcapping.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... D640PA100	Partisil SCX	SCX	10 µm	irregular	85 Å	350 m <sup>2</sup> /g	10	no	Whatman

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L9

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... <b>K140ZX050</b>	Zorbax SCX	SCX	5 µm	spherical	300 Å	50 m <sup>2</sup> /g	-	no	Agilent Technologies

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## 2.1.14 Cation Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### Zorbax SCX

The classical Zorbax product line is typical for family of spherical first generation (type A) packings. The stationary phases show a high hydrophobicity and the support has an excellent mechanical strength. AfC offers the complete product line of the classical Zorbax supports in all column dimensions.

· Bulk material not available.

Zorbax is a registered trademark of DuPont.

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L9

### 2.1.14 Cation Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

#### Sphere-Image 80 SCX

The technical specifications of the Sphere-Image 80 series of packing materials are nearly equivalent to those of Spherisorb 80, a high-performance packing material which has existed for decades and is very well-known. The Sphere-Image 80 packing material is also available with all of the typical modifications. The cation and anion exchangers Sphere-Image SCX and SAX are also available. The most popular materials in the Sphere-Image series are ODS 1 (not endcapped) and endcapped ODS 2. All of these packing materials exhibit excellent batch-to-batch reproducibility.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... D140SI050	Sphere-Image SCX	SCX	5 µm	spherical	80 Å	200 m <sup>2</sup> /g	-	no	Phenomenex
.... D140SI100	Sphere-Image SCX	SCX	10 µm	spherical	80 Å	200 m <sup>2</sup> /g	-	no	Phenomenex

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## 2.1.14 Cation Exchanger silica based

2.1 Packed HPLC Columns and Bulk Material

### PolyEncap WCX

PolyEncap stationary phases are a family of bonded phases where the silica is polymeric encapsulated by acrylic acid ester polymers that are additionally covalent bonded to the silica surface. This type of stationary phases offers different selectivities compared to classical bonded brush type stationary phases. The PolyEncap phases show a slightly higher pH stability in basic mobile phases. Three different bonding chemistries are offered: This are two RP packings with C4 and C18 and one weak cation exchanger (WCX). The C4 also called PolyEncap A has a comparable selectivity to a brush type C8 packing, whereas the corresponding C18 (PolyEncap B) is similar to a brush type C18 packing. The WCX packing is able to separate mono valent and bivalent cations in a single run.

The two RP bondings are available in four particle sizes (3 µm, 5 µm, 7 µm and 10 µm) and two pore sizes (120 Å and 300 Å) whereas the WCX packing is based on a 5 µm, 120 Å silica.

Column Code	Column Dimensions	Packing	Phase	Particle Size	Pore Size	Surface Area	%C	Endcap.	Manufacturer
.... F140PE050	125 x 4.0	PolyEncap WCX	Polymeric encapsulated	5 µm	120 Å	250 m <sup>2</sup> /g	8	no	BISCHOFF
.... F140PE050	250 x 4.0	PolyEncap WCX	Polymeric encapsulated	5 µm	120 Å	250 m <sup>2</sup> /g	8	no	BISCHOFF
.... F140PE050	125 x 4.0 PEEK	PolyEncap WCX	Polymeric encapsulated	5 µm	120 Å	250 m <sup>2</sup> /g	8	no	BISCHOFF
.... F140PE050	250 x 4.0 PEEK	PolyEncap WCX	Polymeric encapsulated	5 µm	120 Å	250 m <sup>2</sup> /g	8	no	BISCHOFF

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- Also available as bulk material in 10 g, 100 g and 1000 g quantities.
- Other column dimensions on request.

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
## 2.1.14 IonEx CAT

### IonEx CAT-S and CAT-125

IonEx CAT-S PBDMA (polybutadiene dimethacrylate) and CAT-125 are a cation exchanger developed by Prof. Schomburg using the polymer coating method. In this phases a spherical silica gel is loaded with a copolymer (PBDMA). The result is a weakly acidic cation exchanger which is excellent for the separation of alkali and alkaline earth metals.

Through observation it has been noted that when using an organic acid with conductivity detection, detection by direct UV is also possible.

Due to the very good mechanical stability of the silica gel particles, IonEx CAT-125 packing material can be used up to 400 bar, CAT-S up to 500 bar. The chemical stability is also very good. The only limitations known until now are that the following solvents cannot be used: 100% methanol, tetrahydrofuran, dioxane and dimethylsulfoxide.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... <b>P140IE050</b>	IonEx CAT-S	silica gel loaded with PBDMA	5 µm	spherical	110 Å	-	-	-	VDS optilab
.... <b>P140IE070</b>	IonEx CAT-S	silica gel loaded with PBDMA	7 µm	spherical	110 Å	-	-	-	VDS optilab
.... <b>E140IE050</b>	IonEx CAT-125	silica gel loaded with PBDMA	5 µm	spherical	100 Å	-	-	-	VDS optilab

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