

Thermo Scientific Accucore XL HPLC Columns



Thermo Scientific Accucore XL HPLC Columns

Based on Core Enhanced TechnologyTM using 4 µm solid core particles, AccucoreTM XL HPLC columns allow users of conventional HPLC methods to enjoy performance far beyond that of columns packed with 5 µm, 4 µm or even 3 µm fully porous particles. Very high separation efficiencies using standard HPLC instruments and conditions provide increased peak resolution and lower limits of detection. An ultra-stable packed bed results in exceptionally robust columns that demonstrate excellent retention and response reproducibility.



Accucore XL HPLC columns Ultimate core performance for conventional HPLC methods


- 4 µm solid core particles for all users
- Same system, same method, better results
- Robust, fast and easy to use



The Key Components of Core Enhanced Technology

Solid Core Particles

4 µm diameter particles with a solid core generate very high efficiencies with conventional HPLC methods



Automated Packing Process

Enhanced automated procedures ensure that all columns are packed with the highest quality

Tight Control of Particle Diameter

Enhanced selection process keeps particle size distribution to a minimum and produces high efficiency columns


Advanced Bonding Technology

Optimized phase bonding creates a series of high coverage, robust phases

Features and Benefits of Accucore XL HPLC Columns

- Compatible with conventional HPLC methods
- High resolution
- Sharp, tall peak shape
- Reproducible chromatography
- Long column lifetime


- No need to change methods or invest in new equipment
- Separate difficult to resolve peaks
- Lower limits of detection detect trace levels of analytes
- Confidence in your results



Ultimate Core Performance for Conventional HPLC Methods

Solid Core Advantage

The solid core particle design, tight control of particle diameter and automated packing process used for Accucore XL HPLC columns result in minimized eddy diffusion.



Minimized eddy diffusion results in higher efficiency – sharper taller peaks

Adjusting Conventional HPLC Methods

For users of conventional HPLC methods working in regulated environments there may be regulatory issues to consider when changing columns in order to realise the improvements offered by newer technologies. For example USP (United States Pharmacopeia) General Chapter <621> Chromatography-System Suitability describes the maximum adjustments that can be made to an analysis so that a method still fulfils the requirements of the system suitability test.

Column Parameter	Allowed Change		
Column length	± 70%		
Column internal diameter	± 25%		
Particle size	Reduction of up to 50%; no increase		
Method Parameter	Allowed Change		
Flow rate	± 50%		
Injection volume	System suitability testing (SST) criteria must be met		
Column temperature	± 10%		
Mobile phase pH	±0.2		
UV wavelength	No changes outside manufacturer specifications		
Concentration of salts in buffer	± 10%		
Composition of mobile phase	Minor component adjustment ± 30% or ± 10% absolute, whichever is smaller		

Transferring a method from a column packed with a 5 μ m fully porous material to an Accucore XL 4 μ m HPLC column requires no changes to method parameters and involves only a 25% reduction in particle size – thus meeting the above requirements.

Accucore XL HPLC Columns

Two different phases are available in the Accucore XL HPLC column range, C18 and C8.

Accucore XL phases are characterized using three tests based on the Tanaka testing protocols. This detailed phase characterization allows the retentivity, selectivity and secondary interactions demonstrated by HPLC packing materials under specified conditions to be objectively compared.

T1: Hydrophobic Interactions

	•••••••••••••••••••••••••••••••••••••••		Parameter	Term
O ^m	HR	Hydrophobic Retention	Retention of compounds based on their hydrophobicity	k'
Ø'''	HS	Hydrophobic Selectivity	Separation of compounds that have similar structure, but differ slightly in hydrophobicity	α
des-	SS	Steric Selectivity	Separation of compounds that have similar structure, but differ in shape	α
SiO HX	HBC	Hydrogen Bonding Capacity	Separation related to degree of end capping	α

T2: Secondary Interactions Under Neutral pH


	••••••		Parameter	Term
SiO NH₂X	BA	Base Activity	Peak shape for basic analytes resulting from total silanol activity (all dissociated at pH 7.6)	t _f
MX	С	Chelation	Peak shapes for chelating analytes resulting from silica metal content	t _f
SiO X+ ph7.6	IEX(7.6)	lon Exchange Capacity (pH 7.6)	Separation between basic and neutral compounds resulting from total silanol activity (all dissociated at pH 7.6)	α


T3: Secondary Interactions Under Acidic pH

			Parameter	Term
SiO OOCHX	Al	Acid Interaction	Interactions resulting in poor peak shape for acidic analytes	t _r
SiO + X ph2.7	IEX(2.7)	lon Exchange	Separation between basic and neutral compounds resulting from acidic silanol activity	α

The results of the phase characterizations are shown in the radar plots on the next page.

Accucore XL C18







The carbon loading of Accucore XL C18 provides high retention of non-polar analytes via a predominantly hydrophobic interaction mechanism. The highly retentive nature of the phase means that it can be used to separate a broad range of analytes.

- Optimum retention of non-polar compounds
- Hydrophobic interaction mechanism
- Separates a broad range of analytes

Accucore XL C8



- Similar selectivity to C18 with lower retention
- Recommended for analytes with moderate hydrophobicity



Accucore XL C8 offers lower hydrophobic retention than columns packed with longer alkyl chain length material, such as C18. It is then therefore recommended for analytes with moderate hydrophobicity, or when a less hydrophobic phase provides optimum retention.

Column Formats

Accucore XL HPLC columns are offered in analytical and micro formats. Optimum conditions and ratings are shown in the table below.

Column ID	Optimum Flow Rate	Optimum Injection Volume	Backpressure Rating	Temperature Rating
2.1 mm	0.3 mL/min	2 µL	600 bar	70 °C
3.0 mm	0.6 mL/min	5 µL	600 bar	70 °C
4.6 mm	1.3 mL/min	10 µL	600 bar	70 °C

Analytical and Micro Columns

Accucore HPLC columns are packed into our high pressure hardware. These stainless steel columns are engineered to the highest quality and have a pressure rating of 600 bar.



Guard Cartridges

Guard cartridges are designed to protect your column from particulates introduced from the matrix or instrument and from any strongly retained components in the injected sample.






4 µm Solid Core Particles for all Users


The 4 µm solid core particles used in Accucore XL HPLC columns have been specifically designed to get the optimum chromatographic performance from conventional HPLC instruments.

- Very high efficiencies
- Little decrease in efficiency as flow rate is increased
- Moderate backpressures

Efficiency

Accucore XL HPLC columns generate higher efficiencies than columns packed with 5 μ m and 3 μ m fully porous material – as shown in the van Deemter curve below.






Identical instrument and method conditions for all columns		
150 x 4.6 mm		
50% water : 50% acetonitrile		
30 °C		
ΙμL		
UV at 254 nm (0.1 s rise time, 20 H		
o-xylene		
1 μL JV at 254 nm (0.1 s rise time, 2		

- 75% higher efficiency than 5 µm fully porous
- 50% higher efficiency than 3 µm fully porous

Backpressure

Accucore XL HPLC columns generate reasonable backpressures, moderately higher than fully porous 5 µm and lower than fully porous 3 µm, that are compatible with conventional HPLC instruments.




• Backpressures between those generated by 3 µm and 5 µm fully porous

• Within conventional HPLC instrumentation pressure limit - even at high flow rates

Impedance

Impedance (E) combines retention time, efficiency and backpressure in a single term. Lower impedance values indicate fast and higher efficiency separations performed at lower backpressures.



Accucore XL C1	8, 4 µm
Fully porous C18,	3 µm
— Fully porous C18,	5 µm


Column Dimensions:	150 x 4.6 mm
Mobile Phase:	50% water : 50% acetonitrile
Temperature:	30 °C
Injection Volume:	1 μL
Detection:	UV at 254 nm (0.1 s rise time, 20 Hz)
Sample:	o-xylene

• 78% lower impedance than 5 µm fully porous

• 67% lower impedance than 3 µm fully porous

Resolution


The high chromatographic efficiencies offered by Accucore XL HPLC columns represent tall, narrow peaks. This provides significant advantages in terms of better peak separations (resolution) and lower limits of detection.



27% higher resolution than 5 μm fully porous
11% higher resolution than 3 μm fully porous

Column Dimensions:	150 x 4.6 mm
Mobile Phase A:	water
Mobile Phase B:	acetonitrile
Gradient:	35 to 60% B in 10 min
Flow Rate:	1.0 mL/min
Temperature:	30 °C
Injection Volume:	5 μL
Detection:	UV at 247 nm (0.1 s rise time, 20 Hz)
Sample:	1. uracil (t0) 2. tebuthiuron 3. metoxuron 4. monuron 5. chlorotoluron 6. diuron 7. linuron

Sensitivity



Column	Amount on Column	Average S/N	Limit of Detection (based on S/N = 3)
Fully porous C18, 5 µm	1 ng	4.86	0.62 ng
Fully porous C18, 3 µm	1 ng	5.31	0.56 ng
Accucore XL C18, 4 µm	1 ng	10.61	0.28 ng

Identical instrument and method conditions for all columns			
Column Dimensions:	150 x 4.6 mm		
Mobile Phase A:	water		
Mobile Phase B:	acetonitrile		
Gradient:	35 to 60% B in 7.5 min		
Flow Rate:	1.3 mL/min		
Temperature:	30 °C		
Injection Volume:	1 µL		
Detection:	UV at 247 nm (0.1 s rise time, 20 Hz)		
Sample:	1. uracil (t0) 2. tebuthiuron 3. metoxuron 4. monuron 5. chlorotoluron 6. diuron 7. linuron (each at 1 ng/µL)		

- 120% more sensitive than 5 µm fully porous
- 100% more sensitive than 3 µm fully porous



Same System, Same Method, Better Results

Ibuprofen and Valerophenone (USP)

The following applications show the improvements in performance that Accucore XL HPLC columns offer without any changes in instrument configuration or method conditions.


$\Delta P = 239 \text{ bar}$ Fully porous C18, 5 µm 30 25 20 15 2 mAU N = 10538 plates 10 S/N = 534 5 0 -5 2.5 5.0 7.5 10.0 12.5 17.5 22.5 30.0 15.0 20.0 25.0 27.5 0 Time (min) $\Delta P = 312 \text{ bar}$ Accucore XL C18, 4 µm 30 25 N = 18274 plates 20 S/N = 1202 15 mAU 10 5 0 -5 2.5 5.0 7.5 10.0 12.5 15.0 17.5 20.0 22.5 25.0 27.5 30.0 0.0 Time (min) Columns Accucore XI C184 um 150 x 46 mm Analytes 1. Valeroph

Analytes		Columns:	Accucore XL C 18 4 µm, 150 X 4.6 mm			
1.	Valerophenone		Fully porous C18 5 µm, 150 x 4.6 mm			
2. Ibuprofen		Mobile phase:	66.3:33.7 (v/v) water with phosphoric acid, pH 2.5:methanol			
		Flow rate:	2 mL/min			
		Column temperature:	30 °C			
		Detection:	UV at 214 nm			
		Injection volume:	5 µL			

73% higher efficiency
125% higher sensitivity



Fat Soluble Vitamins



• 105% higher efficiency

Vitamin E acetate

Vitamin K1

5.

6.

• 28% better resolution of critical pair

Column temperature:

Injection volume:

Detection:

• Reduced run time with good resolution at increased flow rate


30 °C

5 μL

UV at 280 nm



Endocrine Disruptors




Injection volume:

5μL



31% better resolution of critical pair
37% narrower peaks
226% higher sensitivity

Triazines




Injection volume:

5 µL

- 54% better resolution of critical pair
- 29% narrower peaks
- 140% higher sensitivity



Catechins



nalytes		Columns:	Accucore XL C8 4 μm, 150 x 4.6 mm		
	Epigallocatechin		Fully porous C	C8 5 μm, 150 x 4.6 mm	
	Catechin	Mobile phase A:	water + 0.1%	formic acid	
Epigallocatechin gallate		Mobile phase B:	methanol + 0.1% formic acid		
	Epicatechin	Gradient:	Time (min)	% B	
	Gallocatechin gallate		0.0	20	
	Epicatechin gallate		15.0	50	
	Catechin gallate		15.1	20	
	0		20.0	20	
		Flow rate:	1 mL/min		
			1.5 mL/min		
		Column temperature:	25 °C		
		Detection:	UV at 280 nm		
		Injection volume:	5 µL		



1. 2. 3. 4. 5. 6. 7.


129% better resolution of critical pair
34% narrower peaks
Reduced run time with good resolution at increased flow rate

Robust, Fast and Easy to Use


Robustness

Accucore XL HPLC columns are extremely robust offering excellent performance over extended use.

Stability-Retention



Stability-Efficiency



Column:	Accucore XL C8 4 µm, 50 x 2.1 mm			
Mobile Phase:	40:60 acetonitrile:20 mM ammonium formate pH3			
Flow Rate:	0.3 mL/min			
Temperature:	30 °C			
Injection Volume:	2 µL			
Detection:	UV at 233 nm			
Sample:	Non-Steroidal Anti Inflammatory Drugs (NSAIDs) ibuprofen, fenoprofen, naproxen			





Stable retention and efficiency over thousands of injections



Productivity

In addition to using established conventional methods, the high efficiencies offered by Accucore XL HPLC columns, across a wide range of flow rates, allow methods to be optimized to reduce run times and increase productivity.





5	6	7	8 9	10			
Minutes							
Calumn		Mathad	t of loot nook	Reduction in time			
Column		Method	t _R of last peak	Reduction in time			
Fully porous	C18, 5 µm	Standard	8.62 min	0%			
Accucore XL	. C18, 4 µm	Standard	6.56 min	24%			
Accucore XL	. C18, 4 µm	Optimized	4.06 min	53%			
Column Dime	ensions:	150 x 4.6 mm	ID				
Mobile Phase A:		water	water				
Mobile Phase B:		acetonitrile	acetonitrile				
		Standard Me	ethod	Optimized Method			
Gradient:		35 to 60% B i	n 10 min	35 to 60% B in 4 mir			
Flow Rate:		1.0 mL/min		1.3 mL/min			
Temperature	:	30 °C					
Injection Vol	ume:	5 µL					
Detection:		UV at 247 nm (0.1 s rise time, 20 Hz)					
Sample:		1. uracil (t0) 2. tebuthiuron 3. metoxuron 4. monuron 5. chlorotoluron 6. diuron 7. linuron					



Run time reduced by over a third with an improvement in performance

Ordering Information

Accucore XL HPLC Columns

Description	Particle Size	Length (mm)	2.1 mm ID	3.0 mm ID	4.6 mm ID
	ore XL C18 4 μm	50	74104-052130	74104-053030	74104-054630
		100	74104-102130	74104-103030	74104-104630
Accucore XL C18		150	74104-152130	74104-153030	74104-154630
		250	74104-252130	74104-253030	74104-254630
Accucore XL C8	4 μm	50	74204-052130	74204-053030	74204-054630
		100	74204-102130	74204-103030	74204-104630
		150	74204-152130	74204-153030	74204-154630
		250	74204-252130	74204-253030	74204-254630

Accucore XL Guard Cartridges (4/pk)

Description	Particle Size	Length (mm)	2.1 mm ID	3.0 mm ID	4.6 mm ID
Accucore XL C18	4 µm	10	74104-012101	74104-013001	74104-014001
Accucore XL C8	4 µm	10	74204-012101	74204-013001	74204-014001

UNIGUARD Direct-Connection Guard Cartridge Holders

Description	2.1 mm ID	3.0 mm ID	4.6 mm ID
UNIGUARD Drop-In Guard Cartridge Holder	852-00	852-00	850-00
Standard Replacement Tip	850-RT	850-RT	850-RT



Resources for Chromatographers

Thermo Scientific Chromatography Columns and Consumables Catalog

This extensive catalog offers 600 pages of proven chromatography tools and product selection guides. Available online, with a robust search tool and optimized for your iPad[®].

Visit www.thermoscientific.com/catalog



Chromatography Resource Center

Our web-based resource center provides technical support, applications, technical tips and literature to help move your separations forward.

Visit www.thermoscientific.com/crc





For more information visit: www.thermoscientific.com/accucoreXL

© 2012 Thermo Fisher Scientific Inc. All rights reserved. iPad is a registered trademark of Apple Inc., registered in the US and other countries. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 332 3331 France +33 (0)1 60 92 48 34 Germany +49 (0) 2423 9431 20 or 21 United Kingdom +44 (0)1928 534110 Japan +81 3 5826 1615

PS20528 E 06/12S

China +86 21 68654588 +86 10 84193588 +86 20 83145199 800 810 5118 India +91 22 6742 9494 +91 27 1766 2352 Australia 1 300 735 292 (free call domestic) New Zealand 0800 933 966 (free call domestic) All Other Enquiries +44 (0) 1928 534 050 Technical Support North America +1 800 332 3331 Outside North America +44 (0) 1928 534 440

