


APPLICATION NOTES

Fast & Precise Bisphenol A (BPA) Method Toxic Substance Found in Consumer Products



Method Conditions

Column: Cogent Bidentate C8™, 4µm, 100Å
Catalog No.: 40008-75P
Dimensions: 4.6 x 75 mm
Mobile phase: A: DI H₂O/ 0.1% formic acid
 B: Acetonitrile/ 0.1% formic acid

Gradient:	time (min.)	%B
	0	30
	2	30
	6	90
	8	90
	9	30

Flow rate: 0.5 mL/min
Peaks: Bisphenol A
Injection Volume: 5 µL
Detection: UV 275 nm
t_r: 0.9 min


Discussion

Bisphenol A (BPA) is a challenging compound for analysis by HPLC. Biological monitoring of this environmental and health toxicant is a necessary process for surveillance as well as risk assessment. As can be seen from the accompanying chromatograms, a Cogent Bidentate C8™ column was an excellent choice for the analysis of BPA. The peak shape is symmetrical with high efficiency.

The repeatability of the analysis is also remarkable as can be seen in the figure, which shows five overlaid injections of the compound. In addition, the method equilibrates rapidly with only 1 minute post time after the gradient.

For more information
www.MTC-USA.com or
technical@hichrom.co.uk

Red Bull® Energy Drink LC-MS Separation of Taurine & Caffeine



Method Conditions

Column: Cogent Diamond Hydride™, 4µm, 100Å
Catalog No.: 70000-15P-2
Dimensions: 2.1 x 150 mm
Solvents: A: DI H₂O / 0.1% formic acid (v/v)
 B: Acetonitrile/ 0.1% formic acid (v/v)

Gradient:	time (min.)	%B	time (min.)	%B
	0	95	7	30
	1	95	8	95
	6	30		

Temperature: 25 °C
Post Time: 3 min
Injection Vol.: 1 µL
Flow Rate: 0.4 mL/min
Detection: ESI – POS - Agilent 6210 MSD TOF mass spectrometer.

Peaks:
 1. Caffeine 195.0877 m/z [M+H]⁺
 2. Taurine 126.0219 m/z [M+H]⁺


t_r: 0.9 min

Discussion

A method combining solid phase extraction with high performance liquid chromatography–electrospray ionization tandem mass spectrometry was developed for the highly sensitive and accurate screening of caffeine and taurine in energy drinks. The matrix pretreatment by SPE sample extraction was essential in the analysis of taurine. The matrix did not interfere with the analysis of caffeine (data not shown). This method has been successfully applied to screening of caffeine and taurine in commercial Red Bull™ samples, and it is valuable to ensure the safety of food and beverages.

For more information
www.MTC-USA.com or
technical@hichrom.co.uk

Anatoxin-a (ANTX-A) Potent Neurotoxin



Method Conditions

Column: Cogent Diamond Hydride™, 4µm, 100Å
Catalog No.: 70000-15P-2
Dimensions: 2.1 x 150 mm
Solvents: A: 50% MeOH/ 50% DI H₂O/ 0.1% formic acid
 B: Acetonitrile/ 0.1% formic acid

Gradient:	time (min.)	%B	time (min.)	%B
	0	70	6	30
	5	30	7	70

Post Time: 5 min
Injection Vol.: 1 µL
Flow Rate: 0.4 mL/min
Temperature: 25 °C
Detection: ESI – POS - Agilent 6210 MSD TOF MS
Sample: **Stock standard solution:** (0.5 mg/mL) of Anatoxin-a fumarate was prepared in DI water and stored at -20 °C.
Working solution: Sample for injection was diluted 1:100 using 50%A/50%B solvent mixture. It was stored in the dark at 4 °C.


Peak: Anatoxin-a, 166.1226 m/z (M + H)⁺
t_r: 0.9 min

Discussion

Anatoxin-a is a low molecular weight, secondary amine and can be difficult to retain with reversed-phase columns. This application note illustrates a very simple method for retention of Anatoxin-a which does not require derivatization. The retention achieved is good, and the excellent repeatability is demonstrated by the overlay of five injections shown in the figure. If desired, a lower retention time of Anatoxin-a can be obtained by changing the starting concentration of solvent B in the gradient to 60%, designing a steeper gradient, or using a shorter column (e.g. 2.1 x 50 mm).

For more information
www.MTC-USA.com or
technical@hichrom.co.uk

Chlormequat & Mepiquat Plant Growth Regulators



Method Conditions

Column: Cogent Diamond Hydride™ 4µm, 100Å.
Catalog No.: 70000-15P-2
Dimensions: 2.1x150mm
Solvents: A: DI water + 20 mM ammonium acetate, pH adjusted to pH 3.3 with formic acid
 B: Acetonitrile

Mobile Phase: 70%A
Post Time: 5 min
Flow Rate: 0.5 mL/min.

Sample Peaks: 1. Chlormequat (CQ) 122.0737 m/z (M)⁺
 2. Mepiquat (MQ) 114.1277 m/z (M)⁺
 Sample stock solutions were purchased from Sigma. Samples for injection were diluted 1:1000 in the mobile phase.

Detection: ESI – pos. - Agilent 6210 MSD TOF mass spectrometer. CQ and MQ are already charged in solution and under ESI conditions the mass spectra show abundant molecular ion (M)⁺.

Discussion


Analysis of CQ and MQ can be challenging since they are extremely hydrophilic and are only weakly retained in reversed-phase HPLC which makes their separation impossible. In addition these compounds are retained too strongly in normal phase (non polar solvents) chromatography. Due to this situation the preferred method of analysis for these compounds was ion-pair chromatography using high ionic strength (> 150 mM) ion-pair reagents. Buffers used are at very high concentration which makes coupling with MS very difficult if not impossible (source contamination). This note shows a new, sensitive and selective LC-MS method with low ionic strength mobile phase for the analysis of CQ and MQ residues. The method can be used in analysis of many samples including food. The selectivity and sensitivity of the method can be increased by using LC-MS-MS instrument and adequate product ions (CQ 122 m/z to 58 m/z and 63 m/z, MQ 114m/z to 98 m/z and 58 m/z). This method has high repeatability (precision) and reliability of the results.

For more information
www.MTC-USA.com or
technical@hichrom.co.uk

APPLICATION NOTES

Separation of Carotenoids

Use of Shape Selectivity to Resolve Similar Compounds



- Peaks:**
1. Astaxanthin
 2. Capsanthin
 3. Lutein
 4. Zeaxanthin
 5. Canthaxanthin
 6. β -Cryptoxanthin
 7. Echinone
 8. 15-cis β -Carotene
 9. 13-cis β -Carotene
 10. α -Carotene
 11. trans β -Carotene
 12. 9-cis β -Carotene
 13. δ -Carotene
 14. Lycopene

Method Conditions

Column: Cogent C30™, 5 μ m, 200Å
Catalog No.: 71030-25P
Dimensions: 4.6 x 250 mm
Mobile Phase: A: 81/15/4 Methanol / MTBE / DI H₂O (v/v)
 B: 6/90/4 Methanol / MTBE / DI H₂O (v/v)

Gradient:

time (min.)	%B
0	0
90	100

Temperature: 20 °C
Flow Rate: 1.0 mL/min
Detection: UV 450 nm
Sample: Reference standards of each analyte in a mixture.


Discussion

Compounds of the carotenoid family are very lipophilic and often have subtle differences in structure. Resolution of these compounds can be difficult with a typical C8 or C18 column because of the similarities. The Cogent C30™ stationary phase on the other hand can further differentiate by analyte shape in addition to reversed phase interactions. At lower temperatures, the long alkyl chains become more rigid and steric effects become significant, leading to greater selectivity.


For more information
www.MTC-USA.com or
technical@hichrom.co.uk

Acrylamide

Easy & Precise Retention



Acrylamide



Overlay of runs from 2 column lots

Method Conditions

Column: Cogent Diamond Hydride™, 4 μ m, 100Å
Catalog No.: 70000-7.5P
Dimensions: 4.6 x 75 mm
Mobile Phase: Acetonitrile/ 0.1% formic acid
Flow Rate: 1.0 mL/min
Injection Vol.: 1 μ L
Sample: 100 mg/L acrylamide in mobile phase diluent.
Peak: 1. Acrylamide
Detection: 205 nm
t₀: 1.0 min


Discussion

Acrylamide is difficult to retain with conventional reversed phase methods due to its polar nature. With the Diamond Hydride™ column however, retention is readily achievable using a simple isocratic mobile phase. The overlay in the figure illustrates lot-to-lot reproducibility of the stationary phase material. Injections using columns of two different lots of material are shown in the figure.

For more information
www.MTC-USA.com or
technical@hichrom.co.uk

Cyanuric Acid and Melamine

Simple ANP Method



Method Conditions

Column: Cogent Diamond Hydride™ 4 μ m, 100Å.
Catalog No.: 70000-15P-2
Dimensions: 2.1 x 150 mm
Solvents: A: DI water + 0.1% acetic acid
 B: Acetonitrile + 0.1% acetic acid

Gradient:

Time	%B
0.0	100
15.00	50


Post Time: 5 min
Flow Rate: 0.4 mL/min. t₀ = 0.8 min
Sample Prep: 1. Cyanuric acid 1.5 microg/mL: NEG (M-H)-128.0102 m/z analyzed in Neg mode, RT = 1.278 min
 2. Melamine 3 microg/mL: POS (M+H)+ 127.0727 m/z analyzed in Pos mode, RT = 11.260 min

Detection: ESI – NEG for cyanuric acid and ESI – POS for melamine: - Agilent 6210 MSD TOF mass spectrometer.
 Note: When an instrument with a dual ion source is used both samples can be analyzed simultaneously.

Discussion


The benefits of using the method in this application note include: 1. Short equilibration time between gradient runs 2. Excellent repeatability of the results. 3. Use of LC-MS detection eliminates the need for derivatization (as is needed with a GC-MS method). 4. Higher flow rates without compromising the efficiency (higher sample throughput possible) 5. Columns used in this method have very long life and are good for over 2000 injections of real life samples. When analyzing melamine alone, a simple isocratic method can be used. When using the Cogent Diamond Hydride™ column melamine is retained at 80% B mobile phase. The resulting peak shape is very symmetrical.

For more information
www.MTC-USA.com or
technical@hichrom.co.uk



Analysis of Urea

Fast, Reproducible HPLC Method



Urea




Figure above show 5 consecutive runs using 100% water mobile phase without any loss of precision.

Method Conditions

Column: Cogent Bidentate C18™, 4 μ m, 100Å
Catalog No.: 40018-15P
Dimensions: 4.6 x 150 mm
Mobile phase: A: DI water (ISOCRATIC Run: 100%A)
Flow rate: 0.5 mL/min.
Peak/Sample: Urea (1 mg/mL in DI water)
Injection Volume: 10 microL
Detection: UV 210 nm

Discussion

A Cogent Bidentate C18™ column was successfully used to retain Urea when a 100% DI water mobile phase was used. Since this hydrophobic column was made on TYPE-C Silica™ DI water can be used as the complete mobile phase without fear of phase "dewetting" and subsequent loss of retention time. The presented method has the advantage of being very simple, without time consuming derivatization or reaction steps when compared to other methods using ordinary C18 columns. Urea can be determined in complex mixtures using this very fast, reproducible (see insert on the chromatogram) method which could also be useful in the analysis of biological samples for diagnostic purposes.

For more information
www.MTC-USA.com or
technical@hichrom.co.uk

Cogent™ TYPE-C™ silica LC phases

Cogent™ TYPE-C™ silica LC phases have the ability to retain polar solutes at high concentrations of organic solvent by aqueous normal-phase (ANP) and non-polar compounds under reversed-phase (RP) conditions. These revolutionary columns use patented bonding technology to create a surface populated by silicon-hydride functional groups instead of silanols. The lack of surface silanols leads to fast equilibration times, excellent peak shape and extended column lifetimes for a wide range of analytes. These application notes demonstrate the unique abilities of Cogent TYPE-C silica LC columns for a range of clinical analysis applications. Further application notes are available at www.MTC-USA.com or from Hichrom Limited at technical@hichrom.co.uk



Cogent TYPE-C columns can be operated in 3 modes of chromatography: reversed-phase (RP), normal-phase (NP) and aqueous normal-phase. The surface silanols that are present in all Type A and B silicas, even after bonding and extensive endcapping, form a strong association with water resulting in a 'hydration shell' surrounding the silica. However, the silica hydride particles of TYPE-C silica are only slightly hydrophobic and therefore have a weak attraction for water allowing them to be used in aqueous normal-phase (ANP) mode, which unlike HILIC, does not require a 'water-rich' environment in order to operate.

Aqueous Normal Phase (ANP) and Reversed-Phase (RP) Separations




Figure 1. Dual RP and ANP retention capability

Cogent TYPE-C silica based phases (Bidentate C18, Bidentate C8, UDC-Cholesterol, Diamond Hydride, Phenyl Hydride, UDA, Diol and Silica-C) have the ability to operate in ANP mode which enables the retention of polar solutes at high concentrations of the organic component whilst maintaining an aqueous component in the eluent. The exact point in the composition of the eluent where ANP retention begins depends on the solute as well as the stationary phase. In addition, TYPE-C columns can also retain non-polar compounds based on a typical reversed-phase mechanism. Figure 1 illustrates the dual retention capability for both polar (metformin) and non-polar (glyburide) compounds. In this case, with an eluent composition of less than 70% acetonitrile, glyburide and metformin are both retained by a reversed-phase mechanism, with the metformin eluting first. With increasing percentages of acetonitrile, the retention of metformin increases significantly due to ANP mechanisms and now elutes after glyburide.

For further technical advice and additional application notes on Cogent TYPE-C Silica LC columns, contact MicroSolv Technologies, USA, www.MTC-USA.com or global distributor Hichrom Limited, UK www.hichrom.co.uk, technical@hichrom.co.uk