



HICHROM

Chromatography Columns and Supplies

LC COLUMNS
Thermo Scientific
Hypersil

Catalogue 9

Hichrom Limited

1 The Markham Centre, Station Road
Theale, Reading, Berks, RG7 4PE, UK

Tel: +44 (0)118 930 3660 Fax: +44 (0)118 932 3484

Email: sales@hichrom.co.uk www.hichrom.co.uk

Hypersil®

- Spherical porous silica
- 3, 5 and 10µm particle sizes

Hypersil® is a traditional 120Å pore size silica, manufactured by Thermo Scientific. The octyl-, cyano- and phenyl-bonded phases are available as both endcapped and non-endcapped materials.

Hypersil Phases

Hypersil Phase	Functional Group	Endcapped	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Load (%)
Silica	-	No	3, 5	120	170	-
SAS (C1)	Methyl	No	5	120	170	2.5
MOS (C8)	Octyl	No	3, 5	120	170	6.5
MOS-2 (C8)	Octyl	Yes	5	120	170	6.5
ODS (C18)	Octadecyl	Yes	3, 5	120	170	10
ODS-2 (C18)	Octadecyl	Yes	3, 5	80	220	11
Phenyl	Phenyl	No	5	120	170	5
Phenyl-2	Phenyl	Yes	5	120	170	5
CPS (CN)	Cyano	No	3, 5	120	170	4
CPS-2 (CN)	Cyano	Yes	5	120	170	4
APS-2 (NH ₂)	Amino	No	3, 5	120	170	2

Ordering Information

Column i.d. (mm)	Column Length (mm)				Guard Cartridges (4/pk)
	50	100	150	250	
Hypersil 3µm					
2.1	xxxxx-052130 £278	xxxxx-102130 £302	xxxxx-152130 £319	xxxxx-252130 £336	xxxxx-012101 ¹ £161
3.0	xxxxx-053030 £278	xxxxx-103030 £302	xxxxx-153030 £319	xxxxx-253030 £336	xxxxx-013001 ² £161
4.0	xxxxx-054030 £266	xxxxx-104030 £288	xxxxx-154030 £304	xxxxx-254030 £319	xxxxx-014001 ³ £161
4.6	xxxxx-054630 £266	xxxxx-104630 £288	xxxxx-154630 £304	xxxxx-254630 £320	xxxxx-014001 ³ £161
Hypersil 5µm					
2.1	xxxxx-052130 £253	xxxxx-102130 £276	xxxxx-152130 £289	xxxxx-252130 £305	xxxxx-012101 ¹ £146
3.0	xxxxx-053030 £253	xxxxx-103030 £276	xxxxx-153030 £289	xxxxx-253030 £305	xxxxx-013001 ² £146
4.0	xxxxx-054030 £241	xxxxx-104030 £261	xxxxx-154030 £277	xxxxx-254030 £289	xxxxx-014001 ³ £146
4.6	xxxxx-054630 £241	xxxxx-104630 £261	xxxxx-154630 £277	xxxxx-254630 £289	xxxxx-014001 ³ £146

¹ 10 x 2.1 mm. Use with Uniguard® direct connect holder 852-00 (£75) ² 10 x 3 mm. Use with Uniguard® direct connect holder 852-00 (£75) ³ 10 x 4 mm. Use with Uniguard® direct connect holder 850-00 (£75)

When ordering, please replace 'xxxxx' with the appropriate Hypersil material code.

	3µm	5µm		3µm	5µm
Hypersil Silica	30003	30005	Hypersil Phenyl	-	30905
Hypersil SAS	-	30505	Hypersil Phenyl-2	-	31905
Hypersil MOS	30203	30205	Hypersil CPS	30803	30805
Hypersil MOS-2	-	30305	Hypersil CPS-2	-	31805
Hypersil ODS	30103	30105	Hypersil APS-2	30703	30705
Hypersil ODS-2	31603	31605			

Please note that not all phases are available in every dimension.
Please contact Hichrom for ordering information on 10µm Hypersil columns.